

**ARG20671**  
**Goat F(ab')<sub>2</sub> IgG Isotype Control antibody (APC)**

Package: 50 tests

Store at: 4°C

### Summary


Product Description	APC-conjugated Goat Polyclonal antibody as a negative control antibody for Goat F(ab') <sub>2</sub> IgG
Tested Application	FACS, FLISA
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	Goat Fab'2 IgG
Conjugation	APC

### Application Instructions

Application table	Application	Dilution
	FACS	10 µl/10 <sup>6</sup> cells
	FLISA	Assay-dependent
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	

### Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG20671 Goat F(ab')<sub>2</sub> IgG Isotype Control antibody (APC) FLISA image

FLISA: The plate was coated with Rabbit anti-Goat IgG Fab antibody and [ARG23851](#) Mouse anti-Goat IgG Fc antibody [SB115h]. Serially diluted ARG20671 Goat F(ab')<sub>2</sub> IgG Isotype Control antibody (APC) was captured and fluorescence intensity quantified.