

## ARG21792 Mouse anti-Chicken IgY antibody [G-1] (PE)

Package: 100 µg  
Store at: 4°C

### Summary

Product Description	PE-conjugated Mouse Monoclonal antibody [G-1] recognizes Chicken IgY
Tested Reactivity	Chk
Tested Application	FACS, FLISA, ICC/IF, IHC-Fr
Host	Mouse
Clonality	Monoclonal
Clone	G-1
Isotype	IgG1, kappa
Target Name	IgY
Species	Chicken
Immunogen	Affinity purified chicken Ig or isolated lymphocytes
Conjugation	PE


### Application Instructions

Application table	Application	Dilution
	FACS	Assay-dependent
	FLISA	< 1 ug/ml
	ICC/IF	Assay-dependent
	IHC-Fr	Assay-dependent

**Application Note** \* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

### Properties

Form	Liquid
Buffer	PBS, 0.1% Sodium azide and Sucrose.
Preservative	0.1% Sodium azide
Stabilizer	Sucrose
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG21792 Mouse anti-Chicken IgY antibody [G-1] (PE) FLISA image

FLISA: The plate was coated with purified chicken IgY and IgM. Immunoglobulins were detected with serially diluted ARG21792 Mouse anti-Chicken IgY antibody [G-1] (PE).