

ARG21948 Donkey anti-Goat IgG (H+L) antibody (Biotin), pre-adsorbed

Package: 500 µg
Store at: 4°C

Summary

Product Description	Biotin-conjugated Donkey Polyclonal antibody recognizes Goat IgG (H+L)
Tested Reactivity	Goat
Tested Application	ELISA, FACS, FLISA, IHC-Fr, WB
Specificity	Reacts with the heavy and light chains of goat IgG
Host	Donkey
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (H+L)
Species	Goat
Conjugation	Biotin

Application Instructions

Pre Adsorbed Mouse and Rat serum proteins.

Application table


Application	Dilution
ELISA	1:5000 - 1:20000
FACS	Assay-dependent
FLISA	Assay-dependent
IHC-Fr	Assay-dependent
WB	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.5 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Images



ARG21948 Donkey anti-Goat IgG (H+L) antibody (Biotin) (pre-adsorbed) IHC-Fr image

Immunohistochemistry: Frozen section of newborn Mouse cartilage tissue stained with [ARG20787](#) anti-Collagen II antibody (pre-adsorbed) (right) followed by ARG21948 Donkey anti-Goat IgG (H+L) antibody (Biotin) (pre-adsorbed), Streptavidin (FITC) and DAPI.