

ARG23785 Goat anti-Mouse IgG (Fab) antibody (HRP), pre-adsorbed

Package: 500 µl
Store at: 4°C

Summary

Product Description	HRP-conjugated Goat Polyclonal antibody recognizes Mouse IgG (Fab)
Tested Reactivity	Ms
Tested Application	ELISA, FLISA
Specificity	Reacts with the Fab region of mouse IgG1, IgG2a, IgG2b, IgG2c, and IgG3.
Host	Goat
Clonality	Polyclonal
Isotype	IgG
Target Name	IgG (Fab)
Species	Mouse
Immunogen	Mouse IgG
Conjugation	HRP

Application Instructions

Pre Adsorbed Mouse IgG Fc, IgM, and IgA.


Application table

Application	Dilution
ELISA	1:4000 - 1:8000
FLISA	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Purification	Affinity purification with immunogen.
Buffer	PBS (pH 7.4) and 50% Glycerol.
Stabilizer	50% Glycerol
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG23785 Goat anti-Mouse IgG (Fab) antibody (HRP) (pre-adsorbed)
ELISA image

ELISA: The plate was coated with purified Mouse IgG Fab and IgG Fc. Immunoglobulins were detected with serially diluted ARG23785 Goat anti-Mouse IgG (Fab) antibody (HRP) (pre-adsorbed).