

ARG23825 Mouse anti-Rat IgG2a antibody [2A8F4] (Cyanine 5)

Package: 50 µg
Store at: 4°C

Summary

Product Description	Cyanine 5-conjugated Mouse Monoclonal antibody [2A8F4] recognizes Rat IgG2a
Tested Reactivity	Rat
Tested Application	ELISA, ELISPOT, FACS, FLISA, IHC-Fr
Host	Mouse
Clonality	Monoclonal
Clone	2A8F4
Isotype	IgG1, kappa
Target Name	IgG2a
Species	Rat
Immunogen	Rat IgG2a hybridoma
Conjugation	Cyanine 5

Application Instructions

Application table	Application	Dilution
	ELISA	Assay-dependent
	ELISPOT	Assay-dependent
	FACS	< 0.1 µg/10 ⁶ cells
	FLISA	< 1 µg/ml
	IHC-Fr	Assay-dependent

Application Note * The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.

Properties

Form	Liquid
Buffer	PBS and 0.1% Sodium azide.
Preservative	0.1% Sodium azide
Concentration	0.1 mg/ml
Storage instruction	Aliquot and store in the dark at 2-8°C. Keep protected from prolonged exposure to light. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.