

ARG30279
IGFBP4 ELISA Antibody DuoPackage: 1 pair
Store at: -20°C**Component**

Cat. No.	Component Name	Host clonality	Reactivity	Application	Package
ARG10611	anti-IGFBP4 antibody [IBP3]	Mouse mAb	Hu	ELISA	50 µg
ARG10612	anti-IGFBP4 antibody [IBP163]	Mouse mAb	Hu	ELISA	50 µg

Summary

Target Name	IGFBP4
Alternate Names	IGFBP4 ELISA antibody; IGFBP4 antibody

Properties

Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C or below. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Full Name	ELISA Antibody Duo for IGFBP4
Highlight	Related Product: anti-IGFBP4 antibody;