

## ARG40111 anti-PHF11 antibody

Package: 100 µl  
Store at: -20°C

### Summary

Product Description	Rabbit Polyclonal antibody recognizes PHF11
Tested Reactivity	Hu, Ms, Rat
Tested Application	ICC/IF, IHC-P, WB
Host	Rabbit
Clonality	Polyclonal
Isotype	IgG
Target Name	PHF11
Species	Human
Immunogen	Recombinant fusion protein corresponding to aa. 1-292 of Human PHF11 (NP_001035533.1).
Conjugation	Un-conjugated
Alternate Names	NY-REN-34; APY; IGHF; IGER; Renal carcinoma antigen NY-REN-34; NYREN34; PHD finger protein 11; IGHF; BCAP; BRCA1 C-terminus-associated protein

### Application Instructions

Application table	Application	Dilution
	ICC/IF	1:50 - 1:200
	IHC-P	1:50 - 1:200
	WB	1:500 - 1:2000
Application Note	* The dilutions indicate recommended starting dilutions and the optimal dilutions or concentrations should be determined by the scientist.	
Positive Control	SW480	
Observed Size	36 kDa	


### Properties

Form	Liquid
Purification	Affinity purified.
Buffer	PBS (pH 7.3), 0.02% Sodium azide and 50% Glycerol.
Preservative	0.02% Sodium azide
Stabilizer	50% Glycerol
Storage instruction	For continuous use, store undiluted antibody at 2-8°C for up to a week. For long-term storage, aliquot and store at -20°C. Storage in frost free freezers is not recommended. Avoid repeated freeze/thaw cycles. Suggest spin the vial prior to opening. The antibody solution should be gently mixed before use.

Bioinformation


Gene Symbol	PHF11
Gene Full Name	PHD finger protein 11
Function	Positive regulator of Th1-type cytokine gene expression. [UniProt]
Calculated Mw	38 kDa
Cellular Localization	Nucleus. [UniProt]

Images




ARG40111 anti-PHF11 antibody ICC/IF image

Immunofluorescence: A549 cells stained with ARG40111 anti-PHF11 antibody.



ARG40111 anti-PHF11 antibody IHC-P image

Immunohistochemistry: Paraffin-embedded Mouse lung stained with ARG40111 anti-PHF11 antibody at 1:100 dilution.



ARG40111 anti-PHF11 antibody WB image

Western blot: 25 µg of SW480 cell lysate stained with ARG40111 anti-PHF11 antibody at 1:1000 dilution.