

ARG70360
SARS-CoV-2 Nucleocapsid recombinant protein (Omicron B1.1.529 Variant) (Tag Free)

Package: 100 µg

Store at: -20°C

Summary

Product Description	HEK293 expressed, SARS-CoV-2 nucleocapsid recombinant protein (Omicron B1.1.529 Variant)
Tested Application	SDS-PAGE
Target Name	SARS-CoV-2 Nucleocapsid protein (Omicron)
Species	Virus
A.A. Sequence	Met1 - Ala416
Expression System	HEK293

Properties

Form	Powder
Purity	> 95% (by SDS-PAGE)
Buffer	PBS containing 0.5 M NaCl (pH 7.4)
Reconstitution	