

## ARG80374 Human Cardiolipin antibody IgG / IgM ELISA Kit

Package: 96 wells

Store at: 4°C

### Summary

Product Description	ARG80374 Human Cardiolipin antibody IgG / IgM ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human Cardiolipin Ab IgG / IgM in serum and plasma (EDTA, heparin, citrate).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Cardiolipin antibody IgG / IgM
Conjugation	HRP
Sensitivity	IgG:1 GPI U/ml; IgM:0.5 MPI U/ml
Sample Type	Serum and plasma (EDTA, heparin, citrate).
Standard Range	IgG: 7.5 - 120 GPL - U/ml IgM: 5 - 80 MPL - U/ml
Sample Volume	10 µl

### Application Instructions


Assay Time	30, 15, 15, 5 min (RT)
------------	------------------------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.


### Bioinformation

Highlight	Related products: <a href="#">Cardiolipin antibodies;</a> <a href="#">Cardiolipin ELISA Kits;</a> New ELISA data calculation tool: <a href="#">Simplify the ELISA analysis by GainData</a>
Research Area	Immune System kit




ARG80374 Human Cardioliipin antibody IgG / IgM ELISA Kit ELISA image

Patient samples containing high level of specific antibody were serially diluted in samples buffer to demonstrate the dynamic range of the assay and the upper/lower end of linearity. Activity for each dilution was calculated from the calibration curve using a 4-Parameter-Fit with lin-log coordinates.



ARG80374 Human Cardioliipin antibody IgG / IgM ELISA Kit standard curve image

ARG80374 Human Cardioliipin antibody IgG / IgM ELISA Kit results of a typical standard run with optical density reading at 450 nm.



ARG80374 Human Cardioliipin antibody IgG / IgM ELISA Kit standard curve image

ARG80374 Human Cardioliipin antibody IgG / IgM ELISA Kit results of a typical standard run with optical density reading at 450 nm.