

ARG81410 Human CXCL6 ELISA Kit

Package: 96 wells
Store at: 4°C

Component

Cat. No.	Component Name	Package	Temp
ARG81410-001	Antibody-coated microplate	8 X 12 strips	4°C. Unused strips should be sealed tightly in the air-tight pouch.
ARG81410-002	Standard	2 X 10 ng/vial	4°C
ARG81410-003	Standard/Sample diluent	30 ml (Ready to use)	4°C
ARG81410-004	Antibody conjugate concentrate (100X)	1 vial (100 µl)	4°C
ARG81410-005	Antibody diluent buffer	12 ml (Ready to use)	4°C
ARG81410-006	HRP-Streptavidin concentrate (100X)	1 vial (100 µl)	4°C
ARG81410-007	HRP-Streptavidin diluent buffer	12 ml (Ready to use)	4°C
ARG81410-008	25X Wash buffer	20 ml	4°C
ARG81410-009	TMB substrate	10 ml (Ready to use)	4°C (Protect from light)
ARG81410-010	STOP solution	10 ml (Ready to use)	4°C
ARG81410-011	Plate sealer	4 strips	Room temperature

Summary

Product Description	ARG81410 Human CXCL6 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human CXCL6 in serum, plasma (heparin, EDTA) and cell culture supernatants.
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	There is no detectable cross-reactivity with other relevant proteins.
Target Name	CXCL6
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	15.6 pg/ml
Sample Type	Serum, plasma (heparin, EDTA) and cell culture supernatants.
Standard Range	31.2 - 2000 pg/ml
Sample Volume	100 µl

Precision	Intra-Assay CV: 5.2% Inter-Assay CV: 5.6%
Alternate Names	GCP2; SCYB6; GCP-2; Chemokine alpha 3; Granulocyte chemotactic protein 2; Small-inducible cytokine B6; C-X-C motif chemokine 6; CKA-3

Application Instructions

Assay Time	~ 5 hours
------------	-----------


Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol	CXCL6
Gene Full Name	chemokine (C-X-C motif) ligand 6
Function	Chemotactic for neutrophil granulocytes. Signals through binding and activation of its receptors (CXCR1 and CXCR2). In addition to its chemotactic and angiogenic properties, it has strong antibacterial activity against Gram-positive and Gram-negative bacteria (90-fold-higher when compared to CXCL5 and CXCL7). [UniProt]
Highlight	Related products: CXCL antibodies ; CXCL ELISA Kits ; CXCL Duos / Panels ; CXCL recombinant proteins ; New ELISA data calculation tool: Simplify the ELISA analysis by GainData

Images



ARG81410 Human CXCL6 ELISA Kit standard curve image

ARG81410 Human CXCL6 ELISA Kit results of a typical standard run with optical density reading at 450 nm.