

**ARG82062**  
**Mouse CXCL15 ELISA Kit**Package: 96 wells  
Store at: 4°C

### Summary

Product Description	Mouse CXCL15 ELISA Kit is an Enzyme Immunoassay kit for the quantification of Mouse CXCL15 in serum, plasma and cell culture supernatants.
Tested Reactivity	Ms
Tested Application	ELISA
Target Name	CXCL15
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	62.5 pg/ml
Sample Type	Serum, plasma and cell culture supernatants.
Standard Range	125 - 8000 pg/ml
Sample Volume	100 µl
Alternate Names	Scyb15; weche; Small-inducible cytokine B15; C-X-C motif chemokine 15; Lungkine; lungkine

### Application Instructions


Assay Time	4.5 hours
------------	-----------

### Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

### Bioinformation

Gene Symbol	Cxcl15
Gene Full Name	chemokine (C-X-C motif) ligand 15
Function	Chemotactic for neutrophils. Involved in lung-specific neutrophil trafficking during normal and inflammatory conditions. [UniProt]
Highlight	Related products: <a href="#">CXCL antibodies</a> ; <a href="#">CXCL ELISA Kits</a> ; <a href="#">CXCL Duos / Panels</a> ; <a href="#">CXCL recombinant proteins</a> ; New ELISA data calculation tool: <a href="#">Simplify the ELISA analysis by GainData</a>



ARG82062 Mouse CXCL15 ELISA Kit standard curve image

ARG82062 Mouse CXCL15 ELISA Kit results of a typical standard run with optical density reading at 450 nm.