

ARG82147 Acetaldehyde Assay Kit (Colorimetric)

Package: 100 tests
Store at: -20°C

Summary


Product Description	ARG82147 Acetaldehyde Assay Kit (Colorimetric) is a detection kit for the quantification of Acetaldehyde in serum, plasma, urine, tissue, cell culture supernatants, food and beverage.
Tested Reactivity	Hu, Ms, Rat, All
Tested Application	FuncSt
Target Name	Acetaldehyde
Conjugation Note	Read at 565 nm
Sensitivity	2 µM
Detection Range	4 µM - 4 mM
Sample Type	serum, plasma, urine, tissue, cell culture supernatants, food and beverage
Sample Volume	10 µl

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	30 min

Properties

Form	Liquid
Storage instruction	Store the kit at -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.



ARG82147 Acetaldehyde Assay Kit (Colorimetric) typical data demonstration image

ARG82147 Acetaldehyde Assay Kit (Colorimetric) results of a typical data with optical density reading at 565 nm.