

ARG82263 Cellular Senescence Staining Kit (Fluorometric)

Package: 10 assay
Store at: 4°C

Summary

Product Description	ARG82263 Cellular Senescence Staining Kit (Fluorometric) can be used to measure activity of senescence-associated beta Galactosidase in cells or tissue samples by fluorescence microscopy or flow cytometry.
Tested Reactivity	All
Tested Application	FACS, FuncSt, ICC/IF
Target Name	senescence-associated beta Galactosidase
Conjugation Note	λ ex/em = 485 / 520 nm (fluorometric)
Sample Type	Cells or tissue samples.
Sample Volume	Each kit provides sufficient quantities to perform up to 10 assays in 35 mm dishes.

Application Instructions

Application Note	Please note that this kit does not include a microplate.
Assay Time	6 hours to overnight.

Properties

Form	Liquid
Storage instruction	Store the kit at 2-8°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.