

ARG82310 Human Strongyloides IgG / IgM ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG82310 Human Strongyloides IgG / IgM ELISA Kit is an Enzyme Immunoassay kit for the determination of Human Strongyloides IgG / IgM in serum and plasma (citrate, heparin).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Strongyloides IgG / IgM
Conjugation	HRP
Conjugation Note	Substrate: TMB and read at 450 nm.
Sensitivity	Diagnostic sensitivity: 89.47%
Sample Type	Serum and plasma (citrate, heparin). (1:101 diluted)
Standard Range	Cut - off
Sample Volume	100 µl
Precision	Intra-Assay CV: 5.2% Inter-Assay CV: 10.1%

Application Instructions

Assay Time	~ 2 hours
------------	-----------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.