

ARG82648 Human Complement Factor B ELISA Kit


Package: 96 wells

Store at: 4°C

Summary

Product Description	ARG82648 Human Complement Factor B ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human Complement Factor B in serum, plasma, cell culture supernatants, saliva, milk, cell lysate and CSF.
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	

Images



ARG82648 Human Complement Factor B ELISA Kit standard curve image

ARG82648 Human Complement Factor B ELISA Kit results of a typical standard run with optical density reading at 450 nm.