


ARG82650 Human C1q / Complement 1q ELISA Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG82650 Human C1q / Complement 1q ELISA Kit is an Enzyme Immunoassay kit for the quantification of Human C1q / Complement 1q in serum, plasma, cell culture supernatants, urine, saliva, milk and CSF.
Tested Reactivity	Hu
Tested Application	ELISA
Specificity	

Images



ARG82650 Human C1q / Complement 1q ELISA Kit standard curve image

ARG82650 Human C1q / Complement 1q ELISA Kit results of a typical standard run with optical density reading at 450 nm.