

ARG82886 Human Measles IgM antibody ELISA Kit

Package: 96 wells
Store at: 4°C

Summary

Product Description	ARG82886 Human Measles IgM antibody ELISA Kit has been designed for the qualitative determination of specific IgM antibodies against Measles in serum and plasma (Citrate, heparin).
Tested Reactivity	Hu
Tested Application	ELISA
Target Name	Measles IgM antibody
Conjugation	HRP
Sensitivity	Diagnostic Sensitivity: 100 %
Sample Type	Serum and plasma (Citrate, heparin).
Standard Range	Cut - off
Sample Volume	5-10 µl

Application Instructions

Assay Time	1 h, 30, 15 min
------------	-----------------

Properties

Form	96 well
Storage instruction	Store the kit at 2-8°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Research Area	Immune System kit; Microbiology and Infectious Disease kit
---------------	--