

ARG83544

arigoQIK™ Mouse TNF alpha ELISA Development Kit

Package: 1 kit(5 plates), 1 kit
(15 plates)
Store at: 4°C, -20°C

Summary

Product Description

ARG83544 arigoQIK™ Mouse TNF alpha ELISA Development Kit, includes Capture antibody, Detection antibody, Standard, and HRP-Streptavidin Solution. This ELISA Development Kit is designed for the development of sandwich ELISA to measure Mouse TNF alpha in serum, plasma and cell culture supernatants.

For other reagents required for [arigoQIK™ ELISA Development Kit](#), please refer [ARG83524 Integral Reagent Kit \(ELISA Development Kit\)](#)

[More about arigoQIK™:](#)

- Optimized capture and detection antibody pairs
- Reduced incubation time and wash cycles
- 2-hour quicker than conventional ELISA process
- 5- and 15-plate packages available

Tested Reactivity

Ms

Tested Application

ELISA

Target Name

TNF alpha

Conjugation

HRP

Conjugation Note

Substrate: TMB and read at 450 nm.

Sensitivity

7.81 pg/ml

Sample Type

Serum, plasma and cell culture supernatants.

Standard Range

15.63 - 1000pg/ml

Sample Volume

50 µl

Alternate Names

TNF; Tumor Necrosis Factor; TNF-Alpha; TNFSF2; TNFA; DIF 2; Tumor Necrosis Factor Ligand Superfamily Member 2; TNF-A; Tumor Necrosis Factor (TNF Superfamily, Member 2); Tumor Necrosis Factor Ligand 1F; Tumor Necrosis Factor-Alpha; Tumor Necrotic Factor Alpha; TNF Superfamily, Member 2; TNF, Macrophage-Derived; TNF, Monocyte-Derived; APC1 Protein; Cachectin; TNLG1F

Properties

Storage instruction

Store components at 4°C or -20°C. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.

Note

For laboratory research only, not for drug, diagnostic or other use.

Bioinformation

Gene Symbol

TNF

Gene Full Name

Tumor Necrosis Factor

Background

This gene encodes a multifunctional proinflammatory cytokine that belongs to the tumor necrosis factor (TNF) superfamily. This cytokine is mainly secreted by macrophages. It can bind to, and thus functions through its receptors TNFRSF1A/TNFR1 and TNFRSF1B/TNFR. This cytokine is involved in the regulation of a wide spectrum of biological processes including cell proliferation, differentiation,

apoptosis, lipid metabolism, and coagulation. This cytokine has been implicated in a variety of diseases, including autoimmune diseases, insulin resistance, psoriasis, rheumatoid arthritis ankylosing spondylitis, tuberculosis, autosomal dominant polycystic kidney disease, and cancer. Mutations in this gene affect susceptibility to cerebral malaria, septic shock, and Alzheimer disease. Knockout studies in mice also suggested the neuroprotective function of this cytokine. [provided by RefSeq, Aug 2020]

Function Key mediator of cell death in the anticancer action of BCG-stimulated neutrophils in combination with DIABLO/SMAC mimetic in the RT4v6 bladder cancer cell line.


Promotes osteoclastogenesis and therefore mediates bone resorption. [UniProt]

Highlight Related news:
[arigoQIK, DIY your sandwich ELISA kits;](#)

PTM Disulfide bond, Glycoprotein, Lipoprotein, Myristate, Phosphoprotein. [UniProt]

Cellular Localization Cell membrane, Membrane, Secreted. [UniProt]

Images



ARG83544 arigoQIK Mouse TNF alpha ELISA Development Kit standard curve image

ARG83544 arigoQIK Mouse TNF alpha ELISA Development Kit results of a typical standard run with optical density reading at 450 nm.
