

ARG83601 Phospholipid Assay Kit

Package: 96 wells
Store at: 4°C, -20°C

Summary

Product Description	ARG83601 Phospholipid Assay Kit is a detection kit for the quantification of Phospholipid.
Tested Reactivity	Other
Tested Application	FuncSt
Target Name	Phospholipid
Conjugation Note	Read at 540 nm.
Sensitivity	0.2 µM
Sample Type	Serum, plasma and tissues or cell lysates.
Standard Range	0.39 - 25 µM
Sample Volume	50 µL

Properties

Form	96 well
Storage instruction	Store the kit at 4°C, -20°C. Keep microplate wells sealed in a dry bag with desiccants. Do not expose test reagents to heat, sun or strong light during storage and usage. Please refer to the product user manual for detail temperatures of the components.
Note	For laboratory research only, not for drug, diagnostic or other use.